

Solubilization and Purification of Active GST-Fusion Proteins

modified from:

Frangioni J.V. & B.G. Neel (1993) Solubilization and purification of enzymatically active glutathione S-transferase (pGEX) fusion proteins.

Analytical Biochemistry. 210: 179-187.

1. Induce your GST protein of choice and pellet the bacteria. Store in -80°C until you want to purify.
2. Resuspend bacteria in STE buffer with 1X PMSF or complete protease inhibitors with 0.5 to 1.5% sarkosyl (N-laurylsarcosine). If in an eppendorph tube I usually use 1 mL as my resuspension volume. I also don't tend to go above 1% sarkosyl.
3. Lyse bacteria by sonication 3 times for 10 seconds at power setting 60 on sonicator and let sit 5 min on ice between each sonication. If a large volume is being lysed (i.e. 15 mL) lysis can be done using a French Press.
4. Spin out cellular debris by centrifugation at 10,000Xg for 10 min at 4C.
5. Transfer the supernatant to a new tube.
6. Dilute lysate 10 fold with STE buffer.
7. Then incubate lysate on glutathione agarose batch column in tube for 1 hr at 4C. If from an eppendorph tube pellet I usually use 40 μL of glutathione agarose.
8. Wash 5 times with ice cold STE buffer and let sit for 5 min on ice before pelleting glutathione agarose beads and removing supernatant. If I am using an eppendorph tube I use 1 mL washes.
9. Elute GST-fusion protein off of column with 25 mM glutathione in STE buffer with the pH adjusted to 7.0 with 10 M NaOH. If from eppendorph tube I use 50 μL of the 25 mM glutathione STE buffer and do 4 elutions. I do not pool the elutions until after SDS-PAGE analysis and Bradford protein assays to ensure that I am not diluting my protein stocks by combining the elutants.
10. Add glycerol to 5-10% to eluted proteins and take an aliquot for SDS-Page electrophoresis to confirm presence of purified protein. Bradford assays can also be done to determine the concentration of protein recovered.

Solutions and Notes:

STE buffer:

10 mM Tris-Cl, pH 8.0
150 mM NaCl
1 mM EDTA

Add DTT to 5 mM immediately prior to use.

Notes:

Some people put 100 $\mu\text{g/mL}$ of lysozyme in their STE buffer to aid in lysis of the bacterial cells.

This method may be used to solubilize 6-His proteins, but a different buffer must be used. STE buffer contains EDTA which will strip nickel and cobalt columns. If purifying a 6-His protein use PBS or another non-EDTA containing buffer. Use the same non-EDTA containing buffer for making the different imidazole containing buffers for elution of 6-His proteins from the nickel or cobalt column. Most 6-His proteins elute off at 125 mM imidazole. We use washes of 5 mM, 25 mM, and 50 mM imidazole. Elution of protein is done with 250 mM imidazole containing buffers.