

Protocol for Efficient Ligation and Cloning of DNA Fragments with 2-bp overhangs.

Ranjan and Rajagopal (2010). Efficient ligation and cloning of DNA fragments with 2bp overhangs. Analytical Biochemistry 402, 91-92.

Note: The authors used this method for NdeI-HindIII insert ligation. So one of the overhangs is a standard 4bp overhang.

2bp overhang enzymes from NEB:

AseI, AsiSI, BfaI, BspDI, BstBI, ClaI, CviAII, CviQI, HhaI, HinP1I, HpaII, MseI, MspI, NarI, NdeI, PacI, PvuI,

Method:

1. Digest and purify vector and insert.
2. Place 150 ng of digested vector in an eppendorph tube
3. Add insert so vector:insert ratio is 1:3.
4. Heat DNA mixture to 55 C for 2 minutes and then flash freeze at -20 C for 10 minutes.
5. Add remaining ligation components.

i.e. if ligation is 10 μ L

x μ L vector +insert
2 μ L 10X T4 DNA Ligase Buffer
3 μ L (24% w/v) PEG 6000
Bring up to 9 μ L with molecular grade water

Add 1 μ L of T4 DNA Ligase (2000 Units).

6. Let ligation reaction sit overnight at room temperature.
7. Transform bacteria (use 4 μ L of ligation for chemically competent Cells.

24% PEG 6000 solution:

2.4 g PEG 6000 add molecular grade water to 10 mL and filter sterilize