

Isolation of Genomic DNA

From Sambrook et al., Molecular Cloning Manual Vol 2, p.9.16-9.19.

1. Grind tissue in liquid Nitrogen and then place in approximately 10 volumes of **extraction buffer** . Shake and submerge material and then incubate for 1 hr at 37 C.
2. Add **Protease K** to a final concentration of 100 μ g/mL. Gently mix the enzyme into the solution.
3. Place at 50 C for 3 hours. Swirl solution periodically.
4. Cool solution to room temperature and add an equal volume of **Tris-equilibrated Phenol (pH 8.0)** and gently mix the two phases by turning tube end over end for 10 minutes. Separate the two phases by centrifugation at 5,000xg for 15 minutes at room temperature.
5. With a wide bore pipette transfer the aqueous phase to a new tube and extract the aqueous phase with **Tris-equilibrated Phenol (pH 8.0)** two more times.

(note: Do not disturb the interface when transferring the aqueous phase)

6. Precipitate the genomic DNA by adding 0.2 vol **10M Ammonium Acetate** and 2 volumes of 100% ethanol. Swirl to mix and the DNA will precipitate out. Centrifuge at 5000xg for 5 minutes at room temperature.
7. Wash the DNA twice with 70% ethanol and store the pellet in an open tube at room temperature until the last traces of ethanol have evaporated. **Do not allow the pellet to completely dry because it will be difficult to resuspend.**
8. Add 1 mL of 10mM TE (pH 8.0) and gently resuspend the DNA.

Store genomic DNA at 4 C. Repeated freeze thaws will shear genomic DNA.

Extraction Buffer:

10 mM Tris-Cl (pH 8.0)
0.1 M EDTA (pH 8.0)
20 μ g/mL Rnase A
0.5 % SDS