

Glassmilk Purification Kit

Making Glass Milk:

1. Take Borosilliacate Pasteur pipettes and break in a paper towel using a hammer.
2. Using a mortar and pistle grind the pipettes into a fine whilte powder under a fume hood. (This takes approximately 20 minutes).
3. Put ground glass into a 50 mL falcon tube and add 15 mL of sterile water and invert several times. Let sit for 5 minutes and take off supernatant and place in 15 mL falcon tube.
4. Set precipitated glass aside and let dry for use in making more glassmilk at a later date.
5. Take 15 ml slurry and aliquot 1.2 mL into 1.5 mL eppendorph tubes.
6. Spin down glass in 1.5 centrifuge tubes and wash with 1 mL of sterile water
7. Repeat step 6 four times.
8. After final wash resuspend glassmilk in 250 μ L of sterile water.
9. Store in -20 C until needed. Thaw and resuspend by vortexing. Use 20 μ L of slurry for glassmilk purifications.

Low-melt Gel Extraction of DNA Using Glassmilk:

1. Excise DNA band from low-melt agarose gel using a clean scapel and place into a 1.5 mL eppendorph tube.
2. Add 3 vols of 6M NaI solution to gel slice and place at 50 C for 10 minutes to dissolve low-melt agarose (I usually use 500 μ L per slice).
3. Place at room temp for 4 minutes and add 20 μ L of glassmilk slurry. Let incubate at room temperature for 10 minutes. Invert once or twice during incubation.
4. Spin down for 1 min at max speed and remove supernatant.
5. Add 1 mL of New Wash solution and resuspend glassmilk by pipetting then spin down for 1 min at max speed and remove supernatant.
6. Wash glassmilk 2 more times as above with 1 mL New Wash solution.
7. After last wash. Spin glassmilk down one more time and remove any residual New Wash Solution.
8. Add 7 μ L of sterile water and resuspend glassmilk and then place at 37 C for 2 minutes.
9. Spin down max speed 1 min and place supernatant into a new eppendorph tube (this elute contains DNA).
10. Repeat step 8 and 9 pooling both ellutants into the new eppendorph tube.

Solutions:

6 M NaI binding solution:

NaI	90 grams
Sodium sulfite	1 gram
water	up to 100 mL

Filter through Whatmann filter paper and store in a opaque or brown bottle.

New Wash Solution:

50 % Ethanol
20 mM Tris-HCl (pH 7.6-8.0)
1 mM EDTA (pH 8.0)
100 mM NaCl