

Plasmid Isolation Scheme

Modified from Ahn S.C., Baek B.S., Oh T., Song C.S., & Chatterjee B.(2000)
Rapid mini-scale plasmid isolation technique for DNA
sequencing and restriction mapping. *Biotechniques*. **29**,
466-468.

1. Inoculate 2mL LB medium (antibiotic) with a bacterial clone, culture overnight with vigorous shaking at 37°C and harvest bacteria by centrifugation at 11000x g for 1 min.
2. Aspirate culture supernatant and resuspend the pellet completely in 100 µL **resuspension buffer**.
3. Add 100 µL **lysis buffer** and mix gently at room temperature.
4. Add 120 µL **neutralization buffer** and mix gently for 3 min at room temperature.
5. Remove bacterial debris by centrifugation at 11000x g for 1 min and transfer supernatant to a fresh 1.5 mL microcentrifuge tube.
6. Add 200 µL isopropanol to precipitate plasmid DNA and mix thoroughly for 1 min at room temperature.
7. Collect DNA pellet by centrifugation at 11000x g for 30 s and aspirate supernatant.
8. Wash DNA pellet with 500 µL 70% ethanol.
9. Collect the DNA pellet by centrifugation at 11000x g for 30 s and aspirate supernatant.
10. Add 100 µL sterile water to dissolve air-dried DNA.

Solutions:

resuspension buffer:

50mM Tris-HCl, pH 8.0
10mM EDTA, pH 8.0
20 µg RNaseA

Lysis buffer

200 mM NaOH
0, 1% SDS

neutralization buffer

3 M potassium acetate, pH 5.5