

## Electrophoresis of RNA in Formaldehyde Gels

### 1. 5X Formaldehyde Gel-Running Buffer:

0.1 M MOPS (pH 7.0)  
40 mM Sodium Acetate  
5 mM EDTA (pH 8.0)

- 20.6 g 3-(N-morpholino)propanesulfonic acid (MOPS)
- 800 mL Diethyl pyrocarbonate (DEPC)  
treated 50 mM sodium acetate
- adjust pH 7.0 with 2 N NaOH
- Add 10 mL DEPC treated 0.5 M EDTA, pH 8.0

**Autoclave: Solution should be light straw color when done.**

### 2. Formaldehyde Gel

For 1% gel (Agar and water will change if % changes all else stays the same)

- 700 mg Agar into 70 mL of DEPC-treated water
- let cool to 60 C
- add 22 mL MOPS running buffer
- add 20 mL 37% formaldehyde  
(112 mL of gel mix ~1% agarose)

### 3. Denaturation of RNA Samples:

RNA (up to 30 µg)	4.5 µL
5 X formaldehyde gel running buffer	2.0 µL
Formaldehyde	3.5 µL
formamide	10.0 µL

**Incubate samples for 15 minutes at 65 C then chill on ice.**

### 4. Add 2 µL of formaldehyde gel loading buffer

50% sterile glycerol  
1 mM EDTA, pH 8.0  
0.25 % bromophenol blue  
0.25 % xylene cyanol FF

**Prerun gel for 15 minutes at 5 V/cm then load samples**